

## Using The Walker

### Step 1

Align the middle of your foot with the back legs of the walker.



### Step 2

Lift the walker up and place it comfortably in front of you, making sure all four legs are on the floor.



### Step 3

Carefully, walk forward one step at a time.



## Specification

Height - Min	75cm (29.5")
Height - Max	93cm (36.5")
Width	64cm (25")
Depth	57cm (22")
Product weight	3kg (6.6lbs)
Max user weight	136kg (21st)

## Warranty Details

Your Drive DeVilbiss brand product is warranted to be free of defects in materials and workmanship for one year from date of purchase.

The device was built to exacting standards and carefully inspected prior to shipment. This warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the customer.

In the event of a defect covered by this warranty, we will at our option supply parts or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, castors and grips, which are subject to normal wear and need periodic replacement.

If you have any queries or require full warranty conditions, please contact your Drive DeVilbiss stockist.

This does not affect your statutory rights.



Drive DeVilbiss Healthcare Ltd  
Whitehall 26 Business Park,  
1 Heathfield Lane, Birkenshaw,  
BD11 2HW

[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

Document code: 10263KDR\_IFU Issued (05/17)

# Travel Walker

With 5" wheels and glides

DeVilbiss  
HEALTHCARE

drive



[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

## Walker Instructions

Easy “pull up” trigger release for opening and closing




### IMPORTANT SAFETY NOTICE:

Your walker features flip up legs that can be attached to the frame of the walker for easy storage. Be sure that when the legs are turned down for use as a walker, the leg pin is inserted into the lock collar as shown above.

## Walker Adjustments

With your arm extended down, handle should be at wrist height.  
With your hand placed on the hand-grip, your elbow should be flexed 20°-30°.



### CAUTION:

Be aware of the following:


- Electrical cords in way
- Liquids on floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door steps
- Gaps in floor

### WARNING:

1. Use wheels on front legs only or on all four legs with appropriate rear wheel brakes.
2. Do not push one side of the walker ahead and then other side without lifting the walker first.
3. Use Drive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
4. Do not use on stairs or escalators.

## Wheelchair to Walker Transfer

Making certain wheels on the wheelchair are locked lift your body by pushing down on the armrests and carefully transfer one hand at a time to the walker.



## wheelchair to walker transfer

(with lower extremity injury to right leg or foot)

Making certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transferring first the left hand, then the right hand to the walker. (If injury to left leg or foot transfer the right hand first, followed by left.)

