

Invacare® | **Banjo**



A Modern Classic

*Invacare **Banjo** incorporates all the best features of a classic, quality four wheel walker with a fresh design and additional adjustment options. It has been built to perform well outdoors and fold easily for storage when not in use.*



Yes, you can.®



Simplicity and Comfort

The anatomical handle grips are comfortable and offer a secure grip for braking, whilst the height of the handle has an improved range from 710 mm to 980 mm. The safety catch has also been upgraded so that the **Banjo** folds easily through a smart and simple locking system.

Multiple Carrying Options

The *Invacare Banjo* has a hard wearing plastic detachable shelf, a crutch and cane holder, and also spacious basket all fitted as standard.

Features



Tray
Made of impact resistant plastic.



Practical basket
Max load: 9kg.



Crutch holder
Keep a cane or crutch on the rollator.

Technical data

				
Banjo	600 mm	690 mm	710 - 980 mm	610 mm
				
Banjo	390 mm	190 mm	9 kg	130 kg
				260 mm

INVACARE® Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Banjo - UK - 11/2010 - 1504396

Colour



Silver

Product Code

Banjo 1452442



Yes, you can.®