It's ability that matters

DASH

Dashi Mg Power Chair

The ultimate in lightweight, folding power chairs.

Compact and airline friendly.





SUPER LIGHTWEIGHT

The magnesium frame results in one of the lightest chairs on the market with a weight of just 17kg, including the battery.



MANUAL FREEWHEEL LEVERS

Manual freewheel levers on each motor enable you to disable the drive system in order to manoeuvre the chair manually.



ATTENDANT CONTROL

The attendant control option allows an attendant or caregiver to control the powerchair with ease.

Dimensions Total width 600mm Total length (with footplate) 1060mm Total length (without footplate) 720mm Total height 900mm Seat width 500mm Seat depth 410mm 520mm Backrest height Front wheels 6" Rear wheels 8"

Other technical data	
Frame	Magnesium
Front wheel	PU Solid
Rear wheel	Pu Solid
Hub	Magnesium
Motor	250w
Battery capacity	10Ah Lithium
Speed	6 km/h (4 mph)
Travel distance	Up to 9 miles
Attendant control	Optional
Spare battery	Optional



Weight infoChair weight 17kg Max user weight 120kg

Key features

- ⊘ Super lightweight magnesium frame



T: 0845 146 0600

E: sales@rhealthcare.co.uk

F: 01246 453881



RHealthcare