P327 XL Seat Lift

Mid-wheel powerchairs



Your new route to independence



P327 XL Seat Lift

P327 XL SL Specification:

Type Class	B/2
Overall length	112cm / 44"
Overall width	69cm/27"
Floor to seat height Lowest to highest elevation position	62-87cm (min-max)
Seat width	52cm / 20.5"
Seat depth	45.7cm / 18"
Maximum carrying capacity	182kgs / 28 stone
Pattery voltage and capacity	12///55/42
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
	,
Maximum speed	6.4km/h / 4mph
Maximum speed Ground clearance	6.4km/h / 4mph 5cm / 2"
Maximum speed Ground clearance Maximum kerb height	6.4km/h / 4mph 5cm / 2" 5cm / 2"
Maximum speed Ground clearance Maximum kerb height Range* - standard battery pack	6.4km/h / 4mph 5cm / 2" 5cm / 2" up to 32km / 20 miles
Maximum speed Ground clearance Maximum kerb height Range*- standard battery pack The maximum safe slope	6.4km/h / 4mph 5cm / 2" 5cm / 2" up to 32km / 20 miles 6°

^{*}Range on full charge and flat ground

Main Features:

- 10 inch scissor lift mechanism
- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Maximum floor to seat height when fully elevated - 87 cm / 34"



Optional Upgrade:

• Larger seat available - 56cm / 22"

Rascal P327 XL Seat Lift Profile:

The Rascal P327 XL Seat Lift is a heavy-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, with durable full suspension & outstanding manoeuvrability. The P327 XL Seat Lift has a range of up to 20 miles and delivers great performance with 55Ah batteries. This model includes a comfortable Captain's seat with a $10^{\prime\prime}$ scissor lift mechanism.

The scissor lift delivers smooth and quiet lift action and reaches a maximum floor to seat height of 87 cm when fully elevated.

 $Local\,Rascal\,Mobility\,Retailer:$







