



Powered Wheelchair with Reclining Backrest and Footplates

FreedomChair Model T3



A cross-brace designed folding powered wheelchair with removable footplates, a reclining backrest and twin batteries with a range of up to 24 km. Renowned FreedomChair build quality.

Features

Seat Width: 18 in

Reclining backrest

Compact cross-brace design

Removable classic footplates

User Weight Limit: 120kg

Large off-road 12.5" rear tyres

Motors: 2 x 200W continuous / 250W peak

Twin Lithium batteries with a range up to 24km

For customers who prefer the access provided by classic removable wheelchair footplates, the T3 offers a powerful twin motor chair that folds very compactly. The natural flex of the cross-brace design keeps all four wheels in contact with the ground for longer periods over uneven terrain. This chair offers several adjustments to suit the user. The height of the foot plates is adjustable, the armrests can be set at two different heights, and the backrest can be tilted to three different angles.

Specifications

Weight: 28kg with 2 batteries

Maximum weight capacity: 120kg

Seat: 18" wide (extendable)

Speed: 4 mph - Offroad: configurable

Length: 103 cm open: 78 cm folded

Width: 62 cm open: 32 cm folded

Height: 95 cm open: 76.5 cm folded

Minimum turning radius: 88 cm

Driving Range: up to 24 km

Batteries: 2x 24V 7.5AH per battery

Wheels: Front: 8" PU tyre: Rear: 12.5" pneumatic

Motors: 2 x 200W continuous / 250W peak
with electromagnetic braking