

NEW



The Next Generation has arrived!

The all-new **Action²_{NG}** has been developed from the tried-and-tested **Action³** concept. The specifications meet the basic needs of customers whilst refusing to compromise on safety or comfort. The next generation **Action²_{NG}** shares the same essential platform as the **Action³**, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action²_{NG}** perfectly fits into the **Action** family.



Yes, you can.®



Features and Options

Light and comfortable

With a weight from 16.2 kilos, the *Invacare Action²_{NG}* fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and being lightweight makes transportation trouble-free. Cross-compatibility with other members of the *Action* family ensures that the *Action²_{NG}* can be easily adapted to meet individual needs.

Quality of life

The tried-and-tested technology of the *Action* family products allow individuals to go wherever they desire. Cross-compatibility with *Action³_{NG}* and *Action⁴_{NG}* make this an ideal, basic, fleet product. Also available in a transit version.



Castor housing

The castor housing is now mounted outside to have more space for the feet when transferring of foot propelling.



Adjustable center of gravity

The newly-designed rear wheel plate can be adjusted by 50 mm in depth and propose 3 positions in height.



Anti-tippers

Provide extra safety on slopes or uneven ground.

Technical data

Self Propulsion	380 / 405 / 430 / 455 / 480 mm	400 / 430 mm	460 / 485 / 510 mm	430 mm	200 - 280 mm (n°1) 230 / 255 mm (n°2)	330 / 460 mm
Transit	380 / 405 / 430 / 455 / 480 mm	400 / 430 mm	485 / 510 mm	430 mm	200 - 280 mm (n°1) 230 / 255 mm (n°2)	330 / 460 mm
Self Propulsion	80° / 0° - 80°	7° / 0° - 30°	SW + 200 mm	920 / 945 / 970 mm	1030 / 1080 mm	820 / 870 mm
Transit	80° / 0° - 80°	7° / 0° - 30°	SW + 145 mm	945 / 970 mm	950 mm	725 mm
Self Propulsion	325 mm	16,2 kg	8,9 kg	125 kg	Crash tested ISO 7176-19*	
Transit	240 mm	14,9 kg	12,1 kg	125 kg		

INVACARE® Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Action2NG - UK - 02/2011 - 1526737

Frame colour



Anthracite Black

Upholstery



Black nylon

All product dimensions are in millimetres to apply EN 12183 and EN 12184 standards.

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.



Please note colors may vary slightly from those displayed above. These photographs are for general representation of product only.



Yes, you can.®