



A lightweight power chair with excellent range



Quick release
Option: Attendant control

The updated Karma 'Falcon' has several exciting new features including; flip back height adjustable armrests, the new VR2 control system and aluminium split rim rear wheels and castors. The split rims allow the tyres to be replaced easily, saving maintenance costs. Infill rear wheels and solid castors save you worrying about punctures. The new height adjustable armrest is contoured for increased comfort.

Another new feature is the optional quick release steering handle which allows the joystick to be moved from front to back in seconds. This is great for users need assistance occasionally and is very easy to control.

The design has been extensively tested to ensure reliability whilst also benefiting from the experience of



▼ Features

- Aluminum-alloy frame
- Quick release batteries
- Foldable backrest
- VR2 PG controller
- 40Ah batteries
- The heaviest component is merely 32 kg
- Kerb climber

▼ Options

- Foldable drink holder
- Attendant controller system - Quick release? ?
- Parallel swing away joystick bracket

Specifications

59.5cm	104cm	40.6cm (16")	F:7.8"x2" R:12.5"x2.25"	180Wx2	40Ahx2	VR2	N/A	95.5cm	90cm	12°	6km	30Km	Aluminum Alloy	115kg 135kg(Pneu.)	58.6kg 32kg

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 6km/hr (3.7mph). Distance, climbing ability and operating time can vary according to condition.

± 1cm ± 0.5kg

Rev.00 AUG.2007 KARMA reserve the right to modify information herein without further notice.