



IM-9088 Self Propel

IM Lightweight Wheelchairs



IM-9089 Transit

i-Lite Plus Wheelchair

The i-Lite Plus is a modern wheelchair which is both stylish and user friendly. The fold down arms help users get close to table and the black wheels and castors look great against the metallic paint.

It folds down quickly and easily for transport and storage.

Specifications

	Self-propel	Transit
Length (Folded)	830mm	780mm
Length (Unfolded)	1035mm	1035mm
Width (Folded)	330mm	330mm
Width (Unfolded)	690mm	690mm
Height (Folded)	800mm	800mm
Height (Unfolded)	970mm	950mm
Net Weight	11.5kg	13.0kg
Gross Weight	13.0kg	14.5kg
Seat Width	18"	18"
Seat Depth	17"	17"
Seat to Ground	21"	21"
Backrest Height	16"	16"
Wheel Size & Type	8" front 24" rear	8" front 12" rear
TyreType	PU Tyre	PU Tyre
Colour	Silver	Silver
Additional Information	Folding backrest	Folding backrest
Max. User Weight	115kg	115kg

