

2GOABILITY PACE II Specification Sheet

2GOability

Height adjustable push handles

Ergonomic handles with loop style hand brakes for easy grip and control

Weighs just 6.4 kg

Adjustable seat height for extra comfort

Crutch holder as standard

Puncture proof tyres



Convenient mesh basket storage to secure essential items



Sleek modern lightweight aluminium frame

Available Frame Colour



Full Specifications

Seat Width	37 cm	Push Handles Width	52 cm
Seat Depth	24 cm	Maximum Basket Weight	5 kg
Seat Height (From Ground)	55 - 65 cm*	Total Weight	6.4 kg
Total Length	75 cm	Maximum User Weight	100 kg (15 Stone)
Total Width	62 cm	Crash Tested	No
Total Height	101 cm		

KEY: * = Adjustable or optional



www.vanosmedical.co.uk