



Excel Yoga

Van Os Medical UK Ltd.

Excel House | Ashbrooke Park | Lincoln Way | Sherburn in Elmet | North Yorkshire | LS25 6PJ

The Excel Yoga

The Excel Yoga is a lightweight fold away 4mph scooter designed with travel and comfort in mind.

This scooter is supplied with a convenient carry case making the Excel Yoga the perfect travelling companion.



Specifications Excel Yoga

Total length:	99 cm	Turning radius:	216 cm (outer radius)
Total width:	53 cm	Ground clearance:	14 cm
Total height:	62 - 102 cm	Max. speed:	4 mph
Folded length:	47 cm	Range:	18 km
Folded width:	46 cm		
Folded height:	61 cm	Batteries:	24V / 11.5 Ah Lithium batteries
		Charger:	5Ah
Seat width:	42 cm		
Seat depth:	40 cm	Weight:	Scooter only: 27,5 kg incl. battery)
Seat height:	40 cm		Carry case only: 5,75 kg
			Removable armrests: 2 kg each
		User weight:	Max. 110 kg / 17,3 stone
		Available colour(s):	 Sand silver
			 Emerald

Foldable

With the Excel Yoga folding into a convenient carry case it is even easier to transport by car. Or better still take on that holiday you have booked.

As no parts need to be removed, the Excel Yoga folds in one swift easy movement ready to be loaded into its carry case for transportation.



Carry case

The Excel Yoga is equipped with a carry case as standard. The carry case has a telescopic handle suitable for negotiating pavements, roads, train station platforms and departure lounges.

It also has a carry handle on the side for loading/unloading into the car, or onto the luggage trolley



Transportable

Thanks to its compact size when folded, the Excel Yoga is the ideal scooter for transportation.



Armrests

Another feature of the Excel Yoga is the armrests. Unlike our competitors portable travel scooters we include them as standard!

The spec: width adjustable, angle adjustable, removable.



Dashboard

Driving a scooter has never been easier. The Excel Yoga's handlebar steering system with its wig-wag style throttle lever means the Yoga can be driven by one hand.

The Excel Yoga has a simple yet effective dashboard. It is comprised of; battery level indicator, power confirmation LED, high/low speed setting, on/off switch.



Height adjustable

The Excel Yoga has a height and angle adjustable tiller section which allows the user to optimise the tiller position for maximum comfort.

The height adjust mechanism has been manufactured with a quick release lever. This type of system enables the user to set the correct position for the hand controls in the quickest time.



Battery pack

The Excel Yoga is supplied with a removable 11.5Ah Lithium battery pack.

The Excel Yoga has an average range of 18km (depending on terrain, weather & user weight).



Off board charging

As an optional extra, the off-board charging station is an excellent addition to the Excel Yoga's charging facility

Remove the battery from your scooter, place it into the off-board charging tray and plug your existing charger.

You can now charge your Excel Yoga battery anywhere!



Excel Yoga



FOR INTERNATIONAL CUSTOMERS

Koperslagerij 3
4651 SK Steenberg (NB)
The Netherlands

Tel: +31 (0) 167 57 30 20
sales@vanosmedical.nl

www.vanosmedical.nl

FOR UK CUSTOMERS

Excel House
Ashbrooke Park
Sherburn in Elmet
North Yorkshire, England, LS25 6 PJ

Phone: +44 (0) 1977 681 400
Sales@vanosmedical.com

www.vanosmedical.com