

# 9TRL *Technical Information*

Junior Sizes - 9TRLJ Adult Sizes - 9TRL

## OVERALL DIMENSIONS - Attendant Propelled

Seat Width Ref.	Chair Open (mm)			Chair Folded (mm)		
	Width	Length	Height	Width	Length	Height
13	535	840	940	290	610	685
15	585	840	940	290	610	685
16	610	840	940	290	610	685
17	635	915	940	290	685	685
18	660	915	940	290	685	685

## 9TRL(J) WHEELBASE RANGE ACROSS ALL MODELS 13" TO 18"

360mm Smallest Junior to 411mm Largest Adult  
Rear Wheel Size and Function

Item	Attendant
Mounting	Fixed
ø mm	315 Only
Tyre	Puncture Free
Operation	Free

## CASTOR INFORMATION 9TRL(J)

Mounting - Fixed M14 & Tyre - MCP Puncture Free  
Castor size  
Diameter 190mm Width 25mm

## CHAIR WEIGHTS INFORMATION (17")

Version 9TRL(J)

Total weight 14.0kg Weight For Lifting 10.5kg

## MAX. OCCUPANT WEIGHT

127kgs

## SEATING AREA DIMENSIONS

Ranger Ref.	Seat Ref. (Inch)		Effective Measurements (mm)	
	Width	Depth	Width	Depth
9TRLJ	13	15	365	415
9TRLJ	15	16	390	440
9TRLJ	16	16	440	440
9TRL	17	17	465	465
9TRL	18	17	490	465

## CORRIDOR WIDTH (mm) FOR TURNING (mid setting)

Model	Though 360° Between Wall
9TLRJ	1090
9TRL	1150

## ISO 7176 IMPACT STRENGTH & DYNAMIC TESTS

Transport compatible, crash tested at 48kmph/30mph  
Approved for use with webbing tie down and occupant restraint system (wtors) not clamps,  
200,000 cycles two drum, 6,666 cycles kerb drop

# AP100 *Technical*

## OVERALL DIMENSIONS - Attendant Propelled

Seat Width Ref.	Chair Open (mm)			Chair Folded (mm)		
	Width	Length	Height	Width	Length	Height
17	635	915	940	290	685	685

## AP100 WHEELBASE

411mm

Rear Wheel Size and Function

Item	Attendant
Mounting	Fixed
ø mm	315 Only
Tyre	Puncture Free
Operation	Free

## CASTOR INFORMATION AP100

Mounting - Fixed M14 & Tyre - MCP Puncture Free  
Castor size  
Diameter 190mm Width 25mm

## CHAIR WEIGHTS INFORMATION (17")

Version AP100)

Total weight 14.0kg Weight For Lifting 10.5kg

## MAX. OCCUPANT WEIGHT

112kgs

## SEATING AREA DIMENSIONS

Ranger Ref.	Seat Ref. (Inch)		Effective Measurements (mm)	
	Width	Depth	Width	Depth
AP100	17	17	465	465

## CORRIDOR WIDTH (mm) FOR TURNING (mid setting)

Model	Though 360° Between Wall
AP100	1150

## ISO 7176 IMPACT STRENGTH & DYNAMIC TESTS

Transport compatible, crash tested at 48kmph/30mph  
Approved for use with webbing tie down and occupant restraint system (wtors) not clamps,  
Passed 200,000 cycles two drum, 6,666 cycles kerb drop tests



As a result of continual product improvement and development, RHealthcare reserve the right to change product specification without prior notice

# 9TRL & AP100

*General purpose folding back attendant propelled chairs*



**RHealthcare**

It's ability that matters

Sheffield Road, Whittington Moor,  
Chesterfield, Derbyshire S41 8NJ  
Tel: 0845 1460600 Fax: 01246 453 881  
sales@rhealthcare.co.uk  
www.rhealthcare.co.uk



www.rhealthcare.co.uk

**RHealthcare**

It's ability that matters

# 9TRL & AP100

*General purpose folding back attendant propelled chairs*

RHealthcare are a major supplier of standard quality wheelchairs to Health Services across the UK and Ireland. Offering high quality functional chairs at competitive pricing complete with a two year warranty.

The TRL range has also established itself as the chair of choice for Care Homes and Institutions due to its durability, longevity and ready availability of parts.

We pride ourselves on our customer focus and after sales service.



## THE RHEALTHCARE 9TRL & AP100 ARE COMPETITIVELY PRICED QUALITY WHEELCHAIRS FOR EVERYDAY USE.

- 2 Year Warranty
- 9TRL, 9TRLJ & AP100- Attendant propelled chair with 12.5" (315mm) rear wheel & puncture free tyres
- 9TRLJ (Junior) available seat widths - 13", 15" & 16"
- 9TRL (Adult) available seat widths - 17" & 18". AP100 17" seat width
- 9TRL Maximum Occupant weight limit 127Kg
- AP100 Maximum Occupant weight limit 112Kg
- 15 degree hinged folding backrest
- Detachable armrests, c/w Tray Fixing Point
- Detachable, swingaway height adjustable footrests
- Moulded footplate with slip resistant tread
- Durable lightweight steel frame
- Standard puncture proof slim line wheels
- Successfully crash tested and meets ISO Standard 7176/19
- CE marked
- Available in Black
- Fast delivery times of volume and single orders (usually ex stock)
- Low Price Bulk Spares available

## RELAX WITH SUPERHEAD

Superhead is the highly adaptive headrest system suitable for all standard wheelchairs.

- Suitable for use as head support in a transportable situation.
- Fully adjustable head support.
- Provides optimum comfort and is easily adjustable.
- Unique folding mechanism allows system to fold with the wheelchair.
- Fits wheelchairs from 36cm - 50cm.
- Durable yet lightweight aluminium construction.
- High quality, fire retardant fabrics for comfort and safety.
- Available in two sizes.

### Product Suitability

Fits all Standard Wheelchairs

- Medium (seat width 36-44cm)
- Large (seat width 42-50cm)

