

ROMA SENA



ROMA *Global Independence*

SHOPRIDER[®] *The World's Leading Scooter Manufacturer*

ROMA SENA



Special Features

- Captain button padded seat and headrest
- Angle adjustable backrest (fold down)
- Height and angle adjustable seat base
- Articulating front and rear wheel axles, help achieve a more comfortable ride
- Designed in easy to dismantle sections (ideal for transporting or storing)
- Two motor mid-wheel drive system enabling a 360 degree turn within it's own length
- Extendable VSI control unit (for right or left hand use)
- Height and width adjustable arms
- Flip up padded armrests
- Direct access to battery and free wheel release mechanism
- Safety belt
- Puncture proof front and rear swivel castors
- Puncture proof mid-wheel tyres
- Flip up, extendable and height adjustable footrest
- Off board charger

Optional Extras (Subject to suitability)

- Basket (rear)
- Crutch and walking stick holder
- Oxygen bottle holder
- Individual footrest



Specifications

Model	S888WNLS	Maximum Climbing Angle	6° / 10%
Class	1/A	Maximum User Weight	100kg / 16st
Overall Length	940mm / 37"	Battery Size 12v x 2	21Ah
Overall Width	570mm / 22.5"	Range*	24km / 15m
Overall Height	940mm / 37"	Weight (with Battery Pack, Seat, Basket)	63kgs / 139lbs
Seat Width	450mm / 17.5"	Maximum Speed	6 km/h / 3.8mph
Seat Depth	440mm / 17"	Wheel Size (Mid)	200mm / 8"
Seat Height (to ground)	490mm - 590mm / 19" - 23"	(Front)	160mm / 6"
Armrest Height	150mm - 225mm / 6" - 9"	(Rear)	100mm / 4"
Backrest Height from Seat	360mm / 14"	Turning Radius	740mm / 29"
Seat to Footrest/Floorpan Height	350mm - 500mm / 14" - 20"	Battery Charger	3amp off board
Height with Back folded Down	650mm / 25.5"	Ground Clearance**	70mm / 3"

* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel. ** Measurement to underside of chassis.

ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB
Telephone : 01656 674488 Facsimile : 01656 674499
Email : sales@romamedical.co.uk

www.romamedical.co.uk

Dealer Stamp