



Excel *G-Explorer*

Van Os Medical UK Ltd.

Excel House | Ashbrooke Park | Lincoln Way | Sherburn in Elmet | North Yorkshire | LS25 6PJ

The Excel G-Explorer

Van Os Medical's greatest source of pride was the introduction of the markets first modern and inclusive lightweight manual wheelchair the G-Lite Pro.

Move forward 10 years and we bring to you the latest incarnation in this series; The G-Explorer.

We have added pneumatic mountain bike tyres which means the potential for terrains you can negotiate has just increased. That coupled with the textured paint finish gives the chair rugged, outdoor look.



Specifications Excel G-Explorer

Total length: 108 cm (incl. footrests)

Total width: Seat width + 22 cm

Total height: 95-97 cm

Seat width: 40, **45** and 50 cm

Seat depth: 42 cm

Seat height: 54 cm

Back height: 42 cm

Total weight: 11,5 kg

User weight: Max. 113 kg

Available colour(s): Silk Silver and Trail Black

Frame choices

The G-Explorer has two frame colour choices. We have continued the success of the G-Lite PRO by using the silk silver, we have added a textured black finish to give more choice to the user.

Available frame colours:

Trail Black

Silk Silver

Available upholstery colours:

Black nylon upholstery with white stitching

Colour options



Foldable & transportable

The G-Explorer is designed so it can be easily folded and disassembled for transportation, meaning it can fit into the boot of almost any car, into narrow train passage ways and even on plane journeys.

A quick release mechanism engaged by pulling up the seat upholstery folds the chair up in one easy motion.



Backrest & push handles

The side frame has been given a new 4 degree seat angle, providing extra comfort and an ergonomic posture. The backrest angle of the G-Explorer has been adjusted slightly to give the user a more natural seating position providing superb lumbar support.

We added height adjustable push handles to make the attendants experience more comfortable too.



SPECIFICATIONS

		G-Basic	G-Entry	G-Eco	G-Lite Pro
	Total length	107 cm	105 cm	107 cm	108 cm
	Total width	64 cm	63 cm	60-80 cm	60-70 cm
	Total height	92 cm	96 cm	90 cm	95 cm
	Total weight	19 kg	16 kg	19,9 kg	11,5 kg
	Max. occupant weight	135 kg	130 kg	130 kg	113 kg
Legrests	Legrests, swing-away	✓	2 ways	2 ways	2 ways
	Legrests, removable	✓	✓	✓	✓
	Legrest attachment		Plastic	Plastic	Plastic
	Ergonomic release			✓	
	Height adjustment	Quick release	Allen key	Allen key	Allen key
Footplate	Adjustable footplate				
	Calf strap/heel strap/calf plate	Calf strap	Heel strap	Heel strap	Heel strap
Armrests	Armrests, swing-away	✓	✓	✓	
	Armrests, removable				
	Armrests, adjustable				✓
	Stabilisation cube	✓	✓		
Armpads	Armpads, 2 depths mounted			✓	
	Armpads, 4 depths adjustable				
	Extra wide comfort pads				
	Aluminium release				
Push handles	Height adjustable push handles			✓	
Backrest	Backrest, foldable				✓
	Backrest, adjustable				
Upholstery	Upholstery, nylon	✓	✓	✓	✓
	Upholstery, extra	Storage pocket			Removable
Back angle adjustment	Back angle adjustment	Plastic			
Rear wheels	Rear wheel hoop	Plastic	Aluminium	Aluminium	Aluminium
	Rims		Aluminium	Aluminium	Aluminium

SPECIFICATIONS

		G-Basic	G-Entry	G-Eco	G-Lite Pro
	Quick Release axles		✓	✓	✓
	High pressure air tyres				
	Leakfree PU tyres	✓	✓	✓	✓
Front wheels	Front wheel size	8" x 1¼"	200 x 40mm	7" x 1¼"	7" x 1¼"
Side frame	Seat depth adjustable				
Wheelbase	Wheelbase adjustment			2 positions	2 positions
	Rear axle height adjustment			3 positions	
	Plastic front fork		3 positions		
	Aluminium front fork			✓	
	Ball head angle steplessly adjustable				
	Widened ball head				
	Back angle adjustment				
	Seat height adjustment				
Cross frame	Cross frame	Single	Single	Single	Single
Seat upholstery	Nylon upholstery	Foam layer	Foam layer	Foam layer	Removable
	Two-piece and depth adjustable				
	Padded seat				
	Straps, foam & memoryfoam				
Brakes	Brakes	Steel	Aluminium	Steel	Aluminium
	Steplessly depth adjustable	✓	✓	✓	✓
	Steppers/anti-tip wheels	✓	✓	✓	✓
	Crutch holder				
	Adjustable and removable headrest				
	Nylon seat cushion with anti-slip layer				
	Drum brakes				

	G-Evolution	G-Lightweight	G-Modular	G-Neos	G-Neos+	G-Logic	G-Explorer
	✓	✓	✓	✓	✓	✓	✓
				✓	✓		✓
	✓	✓	✓			✓	
	200 x 40mm	7" x 1¼"	8" x 1¼"	8" x 1¼"	8" x 1¼"	7" x 1¼"	7" x 1¼"
			6 positions	Stepless	Stepless		
	2 positions	2 positions	6 positions	Stepless	Stepless	2 positions	2 positions
	3 positions	3 positions	6 positions	3 positions	3 positions		
	✓	✓	✓	✓	✓		
			✓	✓	✓		
				✓	✓		
	Wheel position	Wheel position	Wheel position	Wheel position	Wheel position		
	Single	Single	Single	Double	Double	Single	Single
	Foam layer	Foam layer	Re-tensionable			Removable	Removable
			✓	✓	✓		
				✓	✓		
	Plastic	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
	✓	✓	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓	✓	✓
			✓		✓		
			✓		✓		
					✓		

Excel G-Explorer



FOR INTERNATIONAL CUSTOMERS

Koperslagerij 3
4651 SK Steenberg (NB)
Nederland

Tel: +31 (0) 167 57 30 20
sales@vanosmedical.nl

www.vanosmedical.nl

FOR UK CUSTOMERS

Excel House
Ashbrooke Park
Sherburn in Elmet
North Yorkshire, England, LS25 6 PJ

Phone: +44 (0) 1977 681 400
sales@vanosmedical.com

www.vanosmedical.com