

MARYVILLE RISER RECLINER

New to the Riser Recliner range

we are delighted to introduce the Maryville, a stylish dual motor chair with a range of fabric colour options, it's another great addition to the range.



The Maryville Riser Recliner features a dual motor mechanism, enabling the legrest and backrest to be operated independently, providing countless recline positions. Its compact design is ideal for those with a petite to standard physique and for rooms where space is limited. It is available in 3 attractive fabric colour options, Graphite Grey, Chocolate Brown or Pearl.



LL505

www.drivedevilbiss.co.uk

drive

DeVilbiss
HEALTHCARE




Features:

- » Dual motor provides a wide range of seating and footrest positions
- » Backrest and footrests can be raised up or down either independently or simultaneously with each other, this is ideal for watching TV with your feet up
- » The Maryville chair has an impressive recline function, going almost flat
- » The Maryville rise recliner will help you gently rise to your feet, which helps to give users confidence when both entering and exiting the chair
- » Excellent mid and total snooze positions for supreme comfort and relaxation
- » Generous rise height to help users out of their seat
- » Quiet motors provide a smooth and discreet operation
- » Padded upholstery for added comfort and support
- » Fully padded chaise for increased comfort
- » Convenient side pocket for easy storage
- » In-built fail safe battery back-up system which enables the chair to be raised and lowered in the event of a power failure*
- » Comes with a simple to operate 4 button handset
- » Available in 3 colour options - Chocolate Brown, Graphite Grey, Pearl
- » 2 year warranty on frame and electrical components
- » Maximum weight capacity 150kg (23st)
- » Protective, washable headrest and armrest covers available

*Batteries not included and require two 9V PP3 batteries

Height 111cm	Width 80cm	Depth 88cm	Product Net Weight 45kg
Seat Height 49cm	Seat Width 51cm	Seat Depth 47cm	Weight Capacity 150kg (23st)



If you would like to view this product please contact your local stockist.



Handy Storage Pocket



Fully Padded Chaise



Comfortable TV Position



Total snooze position



Comfortable Padded Upholstery



User Friendly Hand Control

Available in 3 Fabric Colours



**Graphite Grey
Code: CLR41GG**



**Chocolate Brown
Code: CLR41CHOC**



**Pearl
Code: CLR41PRL**

All the dimensions are for the armchairs in their normal chair position. Tolerance is ± 25 mm. All units should be placed with the back at least 60cm away from the wall to prevent damage when reclining the chair. Your riser recliner armchair complies with The Furniture and Furnishings (Fire) (Safety) Regulations 1988.