

**PRODUCT INSTRUCTION SHEET**

**REFERENCE** Step Stool (600, 600S)

**DESCRIPTION**

Sturdy nylon coated tubular steel frame, fitted with non-slip rubber ferrules and a fluted rubber covered wooden top.

**TECHNICAL SPECIFICATIONS**

Top size: 16" (405mm) x 12" (305mm)  
Height of step: 6" (150mm) or 9" (230mm)  
600: 9"  
600S: 6"



**USAGE REQUIREMENTS**

1. Do not exceed user weight of 140kg (22 stone).
2. Clean top with soap and water, and metalwork with cream cleaner.

