

# Van Os Medical UK

## Catalogue

[www.vanosmedical.co.uk](http://www.vanosmedical.co.uk)

Follow us



# Contents

---

Introduction .....	Page 3
Wheelchairs .....	Page 5
Scooters .....	Page 12
Powerchairs .....	Page 17
Powerpacks .....	Page 23
Paediatric .....	Page 24
2GOability .....	Page 25

# Introduction

Welcome to Van Os Medical, the leading supplier of intelligent mobility solutions. We want to help people to live independently for longer, despite disabilities. Using our innovative and reliable equipment they can again enjoy daily life and participate in society.

To achieve our aims, we work according to the following core values:

## Quality



## Trust



### Quality

Is a broad concept for us. We don't only provide quality in our products, but quality is also reflected in our stock, distribution and service. We can always deliver quickly and offer the right support. In all areas we put quality above all.

### Trust

On the basis of the quality we deliver, we build trust. This we have proven since 1971. Our dealers trust our products, trust our staff and end users develop trust in our dealers as their product does not let them down.

## Innovation



## Creativity



### Innovation

We continue to innovate new products and services on the market all the time to stay up to date with our customers and market trends.

### Creativity

In order to continue to innovate flexibility and empathy are key. New ideas from our customers are always welcome, we have many products where the innovation has come from well meaning customers.







# Wheelchairs



G-Logic 



## Excel G-Logic The Customisable Wheelchair

Get personal with the modern new take on the old G-Lite Pro. With a lightweight aluminium frame available in four frame colours, including Trail Black and four upholstery colours; the Excel G-Logic offers a lot of room for customisation. With up to 16 variations, any user can make this chair their own. Also available in Self-propel or Transit, making the G-Logic a very ideal modern manual wheelchair.

### Features

- Half-folding back
- Puncture proof tyres
- Height adjustable armrests
- 16", 18" and 20" seat widths
- Available in Self-Propel or Transit
- Interchangeable upholstery colours
- Optional elevating legrests available
- Swing-away and quick release legrests
- Anti-tip wheels as standard (Self-Propel)
- Quick release wheels for transportation (Self-Propel)
- Attendant brakes as standard for both Transit and Self-propel
- Available frame colours: Ice White, Silk Silver, Matt Black (Standard) and Trail Black



## Excel Globetraveller

### The Lightweight Wheelchair

The Excel Globetraveller has been designed with simplicity and functionality in mind. The ultra compact lightweight aluminium frame makes it the perfect wheelchair for frequent travel. Designed for occasional use, this wheelchair offers the traditional range of compact features, for example, seatbelt, half-folding backrest, swing-away footrests and a simple to use attendant brake system that can easily be locked into place.

- Height adjustable footplates
- Easy disassembly for transport
- Flexible padded seat and back upholstery
- Quick release removable swing-away footrests
- Available frame colour: Blue/Silver



## Caremart Carrymate

### The Foldable Wheelchair

The Caremart Carrymate is a portable wheelchair that is easy to use with its simple folding frame. It can be easily folded to fit inside its carry bag, making it nice and secure for transportation as well as storage; thanks to its extremely lightweight frame. This exceptional example of engineering and design makes the Carrymate stand out from your usual lightweight wheelchairs. It is robust and simple to use with excellent manoeuvrability.

- Storage bag standard
- Extremely lightweight
- Foldable leg and armrests
- Equipped with parking brake
- Available frame colour: White



## Excel G-Explorer

### The Off-Road Wheelchair

Never be afraid to explore the outdoor with the Excel G-Explorer. With its pneumatic 24 inch mountain bike tyres for all terrain, this wheelchair is designed for any adventure. Although not crash tested, this wheelchair has been designed to easily fold and disassemble for transportation. With its ability to fit into almost any car boot, long drives are no hassle.

- Rugged all terrain tyres
- Height adjustable push handles and armrests
- Quick release fold in one mechanism engaged by pulling the seat upholstery
- Available frame colours: Silver and Trail Black





## Excel G3 Hyperlite

### The Modular Sleek Wheelchair

Enjoy effortless propulsion with the all new Excel G3 Hyperlite wheelchair. Designed to be configured to individual needs, the new easy glide wheels and hubs, gives this wheelchair a feather touch feel to the user. With height adjustable armrests and removable footrests with easy to use grips, this chair is unique.

- Flip backward and removable armrests
- Removable swing-away footrests
- Seat belt, tension adjustable backrest and Anti-tip wheels as standard
- Available frame colour: Silver

COMING 2020



## Excel G-Modular FB

### The Prescriptive Wheelchair

The Excel G-Modular is an easily adjustable prescriptive chair that gives the impression of being a truly bespoke product. All of the adjustable components on this chair are colour coded in orange so the user and or the attendant have the instant know-how, on adjusting the chair to the users specific requirements.

- The professionals choice
- Easy disassembly for transport
- Multi-configurable modular chair
- Depth and height adjustable armrests
- Available frame colour: Silver



## Excel G-Nexx

### The Tilt-in-space Wheelchair

The most innovative tilt-in-space (0 - 30 degrees) wheelchair, with a lightweight frame of which all dimensions are steplessly adjustable. Developed with postural support in mind and easily adjustable with a Allen key, the G-Nexx provides ultimate comfort all standard in the base frame.

- Quick and easy adjustable frame
- Seat width and depth stepless adjustments
- Memory foam padding for extra comfort
- Multi-functional 2 gas spring tilt mechanism (one adjusts the backrest and one tilts the chair)
- Available frame colour: Silver



## Excel G3 Kids

### The Paediatric Wheelchair

A robust but lightweight children's wheelchair, simple to fold and handy to transport in almost every car. As an additional option the chair is also available with comfort legrests. It is also equipped with height adjustable push handles and swing-away armrests as standard.

- 20 cm armrest height
- Swing-away removable footrests
- Depth and angle adjustable aluminium brakes
- Available frame colour: Yellow



## Excel G5 Modular Kids

### The Unique Paediatric Wheelchair

With multiple configurations from a single frame which can be adjusted accordingly as your child grows, the Excel G5 Modular Kids is the all rounder. Standard features include seat tension straps, backrest angle adjustment from 90 - 105 degrees and 7 backrest tension straps.

- Non-slip seat cushion
- Angle adjustable footplates
- 3 front fork height adjustments
- Depth and angle adjustable aluminium brakes
- Available frame colours: Blue and Pink



## Excel G5 Modular Junior

### The Popular Paediatric Wheelchair

A versatile children's wheelchair derived from the Excel G-Modular, the Excel G5 Modular Kids big brother. Standard features include height adjustable push handles, backrest angle adjustment with 7 backrest tension straps and height adjustable removable armrests.

- Crutch holder
- Removable Anti-tip wheels
- 22" lightweight aluminium rear wheels
- Depth and angle adjustable brakes with extended levers
- Available frame colours: Blue and Pink





# Specifications

	Excel Globetraveller	Caremart Carrymate	Excel G-Logic	Excel G-Explorer	Excel G3 Hyperlite	Excel G-Modular FB	Excel G-Nexx
Seat Width	45 cm	42 cm	40, 45, 50 cm	40, 45, 50 cm	40, 45, 50 cm	40, 45, 50 cm	40, 45, 50, 55 cm*
Seat Depth	43 cm	40 cm	43 cm	43 cm	43 cm	40, 42.5, 45, 47.5, 50 cm	46.5 - 62.5 cm
Seat Height	50 cm	48 cm	50 cm	50 cm	47.5, 50, 52.5 cm	45, 47.5, 50, 52.5, 55 cm	40 - 56 cm*
Backrest Height	42 cm	40 cm	44 cm	42 cm	42.5, 45, 47.5 cm	39, 41.5, 44 cm	60 cm
Total Length	90 cm	82 cm	95 - 104 cm*	104 cm	103 cm	101 cm	119 cm
Total Width	54 cm	48 cm	62, 67, 72 cm	62, 67, 72 cm	57, 62, 67 cm	62, 67, 72cm	67 - 85 cm*
Total Height	92 cm	89 cm	95 cm	95 - 97 cm*	95 cm	97 cm	103 cm
Total Weight	12 kg	7 kg	12 - 15 kg**	15 kg	15 - 17 kg**	15 - 19 kg**	39 kg
Maximum User Weight	100 kg (15 Stone)	100 kg (15 Stone)	113 kg (16", 18") 135 kg (20")	113 kg (16", 18") 135 kg (20")	130 kg (16", 18") 150 kg (20")	130 kg (20 Stone)	158 kg (24 Stone)
Transit	Yes	Yes	Yes	No	Yes	Yes	Yes
Self-propel	No	No	Yes	Yes	Yes	Yes	Yes
Backrest Foldable	Yes	Yes	Yes	No	Yes	Yes	No
Back Angle Adjustment	No	No	No	No	No	Yes	Yes
Attendant Brakes	Yes	Yes	Yes	No	No	Optional	Yes
Anti-tip Wheels	No	No	Yes	Yes	Yes	Yes	Yes
Armrests (Swing-away)	No	Yes	No	No	Yes	Yes	Yes
Armrests (Removable)	No	No	No	No	Yes	Yes	Yes
Armrests (Adjustable)	No	No	Yes	Yes	Yes**	Yes	Yes
Extra Wide Comfort Pads	No	No	No	No	No	No	Yes
Height Adjustable Push Handles	No	No	No	Yes	No	No	Yes
Seat Upholstery (Removable)	No	No	Yes	Yes	Yes	Yes	No
Padded Seat	No	No	Yes	Yes	Yes	Yes	Yes
Legrests (Removable)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Calf Strap/Heel Strap/Calf Plate	Calf Strap	No	Heel Strap	Heel Strap	Heel Strap	Heel Strap	Calf Plate
Crash Tested	No	No	Yes	No	Yes	Yes	Yes
Available Colour Options	Blue/Silver	White	White Silver Matt Black Trail Black	Silver Trail Black	Silver	Silver	Silver

## Excel G3 Kids

## Excel G5 Modular Kids

## Excel G5 Modular Junior

30, 35 cm	30, 32.5, 35 cm	35, 37.5, 40 cm
30, 35 cm	30, 32.5, 35 cm	35, 37.5, 40 cm
40.5, 43 cm	40, 42, 44.5, 45.5, 46, 46.5 cm	40, 42.5, 44, 45.5, 46.5 cm
38, 40, 42 cm	40, 42.5, 45 cm	40, 42.5, 45 cm
83 cm	86 cm	90 cm
51, 56 cm	51, 53.5, 56 cm	56, 58.5, 61 cm
80 cm	84 cm	83 cm
13.9 kg	14.4 kg	14.9 kg
75 kg (11 Stone)	80 kg (12 Stone)	82 kg (13 Stone)
No	No	No
Yes	Yes	Yes
No	No	No
No	Yes	Yes
No	No	No
Yes	Yes	Yes
Yes	Yes	No
No	Yes	Yes
No	Yes	Yes
No	No	No
Yes	Yes	Yes
Yes	Yes	Yes
Yes	Yes	Yes
Yes	Yes	Yes
Yes	Yes	Yes
Heel	Heel	Heel
Yes	Yes	Yes
Yellow	Blue Pink	Blue Pink

**KEY:** **Bold** = Standard option \* = Adjustable or optional \*\* = Varies on prescription options







# Scooters



## Excel Galaxy II

### A Market Leading Scooter

Introducing the new Excel Galaxy II scooter, available in 7 colours, including 2 new and exclusive limited edition matte silk finish colours, the Galaxy II is the definition of modern. With a new multifunction LCD display screen, adjustable steering column, swivel chair and maximum speed of 8 mph, this maxi size road scooter is ideal for an active user. With a mile range varying between 25 - 35 miles, you can explore in ultimate comfort and style, handling all terrains including gravel, grass, pavements and bumpy roads.

### Features



- 360 degree rotating chair
- 90 Ah batteries for larger range
- New and improved LCD tiller display screen
- Available colours: Dark Storm\* and Rapid Ocean\* (Customer Limited Exclusives)\*
- Standard colours: Gun Metal Grey, Shiraz Red, Emerald Green, Lightning Blue and Midnight Black.





## Excel Roadster DX8 Deluxe

### The Deluxe Compact Scooter

The Excel Roadster DX8 Deluxe is just as striking as the original DX8, but incorporating a new 4-pole technology motor, bigger batteries, added suspension and larger 12 inch wheels for increased range. This 8 mph scooter comes with a range of features including pneumatic tyres, comfortable adjustable Captain seat swivel chair, angle adjustable tiller and a new modern two-tone shroud design.



### Features

- 12 inch wheels
- Larger batteries (55Ah standard)
- Full LED lighting system
- Comfortable and compact
- Eye catching two-tone shroud design
- Available colours: Gallardo Orange, Benz Silver, Tuscan Red and Kite Blue



## Excel Roadster DX8

### The Compact Road Scooter

The Roadster DX8 is a striking midi sized 8 mph scooter. Comfortable and compact enough to be stored with simplicity, it has the speed but agility to get you around town and beyond. Featuring a comfortable Captain seat swivel chair, pneumatic tyres for extra comfort and a multi position tiller.

- Dashboard battery indicator
- Swivel chair with folding backrest
- Full LED lighting system and design
- Available colours: Benz Silver, Gallardo Orange, Kite Blue and Tuscan Red



# Travelux Tiempo

## The Transportable Scooter

The Travelux Tiempo is an extremely comfortable and lightweight 4 mph travel scooter. Robust and highly manoeuvrable, this scooter is great for both indoor as well as outdoor use. With its all new and striking modern design, the Tiempo is the ideal solution for very active users.

The design of the scooter makes it very simple to break down into 5 pieces, of which each part has a handle to help control lift. Not only does the Tiempo have a new and improved Delta steering wheel for added comfort, the wrap around style minimises strain on the wrists which allows for both left and right handed operation.



### Features

- Puncture proof tyres
- Delta steering wheel
- LED lighting for better visibility and energy efficient longer lasting bulbs
- 4 mph speed scooter and 10 mile range
- Available colours: Benz Silver, Cooper Blue, Gallardo Orange and Tuscan Red



# Excel Yoga

## The Lightweight Foldable Scooter

Lightweight and fold-away, the Excel Yoga is a 4 mph scooter designed with travel and comfort in mind. The Yoga is named for its ability to fold into a single unit which then allows it to fit conveniently into the provided carry case; the perfect travel companion. Created in mind to give the user freedom to take their mobility scooter with them without fear of damaging their product.

- Folds in one
- Effective dashboard
- Carry case with coded lock
- Height and angle adjustable tiller
- Width adjustable and removable armrests
- Available colours: Sand Silver and Tuscan Red



# Specifications

	Travelux Tiempo	Excel Roadster DX8	Excel Roadster DX8 Deluxe	Excel Yoga	Excel Galaxy II
Total Length	102 cm	120 cm	132 cm	99 cm	146 cm
Total Width	51 cm	61 cm	61 cm	42 cm	65 cm
Total Height	90 cm	99 cm	102 cm	62 - 102 cm*	130 cm
Seat Width	41 cm	46 cm	46 cm	42 cm	51 cm
Seat Depth	40 cm	40 cm	40 cm	38 cm	46 cm
Seat Height (From Ground)	50 - 55 cm*	53 cm	53 - 58 cm*	58 cm	65 cm
Turning Radius	115 cm	153 cm	153 cm	99 cm	161 cm
Ground Clearance	8 cm	9 cm	11 cm	14 cm	13 cm
Maximum Speed	4 mph	8 mph	8 mph	4 mph	8 mph
Range	19 km (10 miles)	35 km (22 miles)	43 km (27 miles)	18 km (11 miles)	56 km (35 miles)
Battery	2 x 12V / 12Ah	2 x 12V / 35Ah	2 x 12V / 55Ah	1 x 24V / 11.5Ah	2 x 12V / 90Ah
Charger	24V / 2Ah	24 V / 4Ah	24V / 4Ah	24V / 5Ah	24V / 8Ah
Total Weight (Including Battery)	48 kg	86 kg	92 kg	27.5 kg	127 kg
Maximum User Weight	113 kg (17.7 Stone)	136 kg (21 Stone)	136 kg (21 Stone)	110 kg (17 Stone)	160 kg (25 Stone)
LCD Screen Display	No	No	No	No	Yes
Adjustable Tiller	Yes	Yes	Yes	Yes	Yes
Delta Tiller	Yes	Yes	Yes	No	Yes
Adjustable Armrests	Yes	Yes	Yes	Yes	Yes
Adjustable Seat Height	Yes	Yes	Yes	No	Yes
Angle Adjustable Backrest	No	Yes	Yes	No	Yes
Swivel Chair	Yes	Yes	Yes	No	Yes
LED Lighting	Yes	Yes	Yes	No	Yes
Available Colour Options	Benz Silver Cooper Blue Gallardo Orange Tuscan Red	Benz Silver Gallardo Orange Kite Blue Tuscan Red	Benz Silver Gallardo Orange Kite Blue Tuscan Red	Sand Silver Tuscan Red	Emerald Green Gun Metal Grey Lightning Blue Midnight Black Shiraz Red Dark Storm** Rapid Ocean**

**KEY:** \* = Adjustable or optional \*\* = Limited Edition *Please Note: Mile range varies depending on weather and terrain conditions.*







# Powerchairs

## Excel Airide S-Preme

### The Unlimited Powerchair

The S-preme is the top of the range model in the Airide range. Building on the extensive features of all the other models, the S-preme offers up to an extra 8 electrical adjustment functions, along with the ability to add driving lights and indicators if desired.

### Features

- Anti-tips as standard
- Swing-away left or right hand joystick
- Removable adjustable and swing-away legrests
- Customisable seating and backrest options from our Silver seating range, tailored to meet each users individual requirements
- Available frame colour: Silver



# THE AIRIDE

## Excel Airide X-Tend

### The Electrical Adjustments Powerchair

Building on the popularity of the Airide B-Ace, the X-Tend is the next step up in the Airide range. Being one of the most popular outdoor and indoor chairs, this chair is fully crash tested and also has two electrical adjustments to the seat position for added comfort. The Airide X-Tend allows the user to specially choose what options are required at point of order or during ownership.

- 20 mile range
- 120 cm turning radius
- 4 - 6 mph speed powerchair
- Swing-away left or right hand joystick
- Height and depth adjustable armrests



## Excel Airide B-Ace

### The Comfortable Powerchair

The Excel Airide B-Ace is the base drive only model, offering excellent features, sturdy frame and also comes with as standard heavy duty 90 Ah batteries. The B-Ace gets its name from one of its best features; the chair is supplied with air-pressure suspension as standard that can be set to your own specific requirements. Constructed using a rigid yet lightweight aluminium frame, this Airide is in a league of its own.

- 20 mile range
- Optional headrest
- Air pressure suspension
- 4 mph speed powerchair
- Swing-away left or right hand joystick



## Excel Airide Compact

### The Compact Powerchair

The entry level of Airide powered wheelchair series, the Excel Airide Compact is the entry level of the Airide powered wheelchair series. Offering excellent stability, this powerchair is most definitely one of the best value for money indoor and outdoor powerchairs on the market today.

- 17 mile range
- 4 mph speed powerchair
- Left or right hand joystick
- Tension adjustable backrest
- Swing-away left or right hand joystick
- Width, depth, height adjustable and removable armrests





## Excel Airide Go!

### The Illuminated Powerchair

The base frame of the Airide Go! has the same driving dynamics and configurability as the Airide Compact. However, the Go! is equipped with 14 inch drive wheels and 10 inch wheel as standard, along with LED lights as well as indicators. In addition, the electronics are ready for instalment of an electrical tilting system.

- 17 mile range
- LED lights and indicators
- 4 mph speed powerchair
- 10 tension adjustable backrest
- Swing-away left or right hand joystick



## Travelux Quest

### The Ideal Compact Powerchair

With the small turning radius the Travelux Quest is the perfect product for use in tight spaces, making it the ideal robust powerchair for both indoor and or outdoor use. This chair can also be easily transported as it breaks down into 4 smaller pieces.

- 50 Ah controller
- Under Seat Basket
- 9" PU rear tyres with designer wheels
- Flip backward and width adjustable armrests
- Angle adjustable footplate (0-35 degrees)
- Available frame colour: Red



## Travelux Quest Midwheel

### The Comfortable Powerchair

The Travelux Quest Midwheel powerchair boasts a smoother travel experience with its ability to rotate around its axis. With the help of its independent swivel frame construction this creates an exceptional dynamic ride.

- 360 degree swivel height adjustable and removable seat
- Motor with electromagnetic brake
- Foldable and width adjustable armrests
- Standard Fish-on seat (Captain seat available)\*
- Available frame colour: White



# Travelux Venture

## The Deluxe Powerchair

Ideal for indoor use, the Travelux Venture is a compact manoeuvrable elevating electric wheelchair with a small turning radius of only 25 inches. With a comfortable elevating seat which allows the user to raise and lower the seat height with the touch of a button. This fantastic new feature will help you to access everyday heights.

- Fold-away one piece footplate
- Width and height adjustable armrests
- Standard adjustable Captain seat with adjustable seat height and elevating function
- Available frame colour: White



# Specifications

	Travelux Quest	Travelux Quest Midwheel	Travelux Venture	Airide Compact	Airide Go!	Airide B-Ace
Total Length	84 cm	98 cm	102 cm	114 cm	114 cm	114 cm
Total Width	57 cm	63 cm	56 cm (Fish-on) 64cm (Captain)	57 cm	57 cm	65 cm
Total Height	84 - 86 cm*	84 - 86 cm*	94 - 99 cm*	99 cm	99 cm	103 cm
Seat Width	43 cm	43 cm	45 cm	42.5 cm	42.5 cm	45 - 55 cm*
Seat Depth	43 cm	43 cm	40 cm	50 cm	50 cm	40 - 52.5 cm*
Seat Height (From Ground)	52 - 54 cm*	52 - 54 cm*	53 - 69 cm*	50 cm	54 cm	43 cm
Backrest Height	40 cm	40 cm	48 cm	55 cm	55 cm	50 cm
Turning Radius	82 cm	82 cm	63 cm	120 cm	120 cm	120 cm
Ground Clearance	4 cm	4 cm	4 cm			
Maximum Speed	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph
Range	20 km (12 miles)	20 km (12 miles)	26 km (16 miles)	27 km (17 miles)	27 km (17 miles)	32 km (20 miles)
Battery	2 x 12V / 18Ah	2 x 12V / 18Ah	2 x 12V / 35Ah	2 x 12V / 50Ah	2 x 12V / 50Ah	2 x 12V / 90Ah
Charger	24V / 3Ah	24V / 3Ah	24V / 3Ah	24V / 5Ah	24V / 5Ah	24V / 8Ah
Total Weight (Including Battery)	67 kg	67 kg	71 kg	118 kg	118 kg	120 kg
Maximum User Weight	135 kg (21 Stone)	135 kg (21 Stone)	135 kg (21 Stone)	160 kg (25 Stone)	160 kg (25 Stone)	160 kg (25 Stone)
Adjustable Armrests	Yes	Yes	Yes	Yes	Yes	Yes
Joystick Contoller	Yes	Yes	Yes	Yes	Yes	Yes
Seatbelt	No	No	No	Yes	Yes	Yes
Angle Adjustable Backrest	No	No	No	Yes	Yes	Yes
Headrest	No	No	No	Yes*	Yes*	Yes*
Adjustable Legrests	No	No	No	Yes*	Yes*	Yes
Adjustable Footplate	Yes	Yes	Yes	Yes	Yes	Yes
LED Lighting	No	No	No	No	Yes	No
Rear Wheel Drive	Yes	No	No	Yes	Yes	Yes
Fish-On-Seat Available	Yes	Yes*	Yes*	No	No	No
Captain Seat Available	No	Yes*	Yes*	No	No	No
Silver Seating	No	No	No	No	No	Yes
Available Colour Options	Red	White	White	Silver	Silver	Silver



Airide X-Tend	Airide S-Preme
114 cm	114 cm
65 cm	65 cm
103 cm	103 cm
45 - 55 cm*	45 - 55 cm*
40 - 55 cm*	40 - 55 cm*
43 - 47 cm*	43 - 47 cm*
50 cm	50 cm
120 cm	120 cm
4 - 6.25 mph*	4 - 6.25 mph*
32 km (20 miles)	32 km (20 miles)
2 x 12V / 90Ah	2 x 12V / 90Ah
24V / 8Ah	24V / 8Ah
120 kg	120 kg
160 kg (25 Stone)	160 kg (25 Stone)
Yes	Yes
Yes	Yes
Yes	Yes
Yes	Yes*
Yes*	Yes*
Yes	Yes
Yes	Yes
Optional	Yes
Yes	Yes
No	No
No	No
Yes	Yes*
Silver	Silver

**KEY:** \* = Adjustable or optional. *Please Note:* Mile range varies depending on weather and terrain conditions. Speed also varies depending on prescription options.

# Powerpacks

## Excel Click & Go 'Lite'

### Just Click and Go!

The smallest and lightest in the Click & Go model range, easy to transport as well as quick to install to most foldable wheelchairs.

With its two wheels give good contact with the ground, offering consistent levels of traction to assist with propelling the chair. Ideal for shorter distances with flat terrain. The Click & Go Lite can be easily removed in seconds; returning the chair to a standard wheelchair. This power-pack is very easy to operate and you can also adjust the maximum speed.

### Features

- PU tyres
- Hand control unit
- Range of 8 km (5 miles)
- Forward and backward drive
- Includes mounting kit as standard



## Specifications

Excel Click & Go 'Lite'

Total Length	36 cm
Total Width	25 cm
Total Height	29.5 cm
Suitable Wheelchair Seat Width	40 - 50 cm
Drive Wheels	20 x 5 cm
Maximum Speed	5.5 km/h
Range	8 km (5 miles)*
Battery	1 x 12V / 20Ah
Charger	12V / 3Ah
Total Weight (Including Battery)	13.2 kg
Maximum User Weight (Including Wheelchair Weight)	135 kg (21 Stone)
Hand Controller	Yes
Fault Diagnosis	Yes
Battery Level Indicator	No
Adjustable Push Handles	No
Adjustable Speed	Yes
Forward Drive	Yes
Reverse Drive	Yes
Electromagnetic Brake	Yes
Circuit Breaker	No
Removable Battery Pack	Yes
Available Colour Options	Silver

*Please Note: \*Mile range varies depending on weather and terrain conditions.*

# Paediatric

## Excel Elise Travel Buggy

### The Lightweight Pushchair

The Excel Elise Travel Buggy is a high quality and lightweight aluminium paediatric pushchair. In just one swift move, you can turn the buggy into a very compact transportable package. Thanks to the range of different optional accessories, you can offer your child excellent protection and comfort in all weather conditions.

### Features

- Rear wheels with foot brake
- Aluminium lightweight frame
- Double front castor wheels for extra stability
- Optional accessories available for purchase\*



## Specifications

### Excel Elise Travel Buggy

Seat Width	40 cm
Seat Depth	27 cm
Seat Height Front	60 cm
Seat Height Back	47 cm
Backrest Width	40 cm
Backrest Height	56 cm
Total Height	112 cm
Total Width	55 cm
Total Length	96 cm
Total Weight	6.9 kg
Maximum User Weight	50 kg (8 Stone)
Available Colour Options	Red/Grey Blue/Grey

### Please Note:

*Elise Buggy pictured with all optional accessories, which are not included as standard\**

## 2GOability Pace II The Lightweight Rollator

The 2GOability Pace II is a lightweight aluminium framed rollator, that offers high stability and simplistic usability. Well manufactured with a range of features designed to help make it a more comfortable walker for both indoor and or outdoor use. The risk of trapping fingers in the frame is eliminated with the Pace II's unique easy folding system.

Features include ergonomic handles and loop type hand brakes to slow or stop the walker. The brakes can also be locked to create a stable platform for the user when resting or sitting throughout daily tasks.

### Features

- Crutch holder as standard
- Removable seat and backrest
- Height adjustable push handles
- Mesh basket for securing items
- Available frame colour: Silver



## 2GOability Access The Affordable Transit Wheelchair

The 2GOability Access wheelchair is a comfortable, lightweight aluminium Transit wheelchair, with a really helpful range of high quality features. With its ultra compact frame and lightweight foldable design, the Access wheelchair is perfect for travel.

- Weighs only 9.5 kg
- Puncture proof tyres
- Foldable for easy transport and storage
- Quick and easy set up (no tools required)
- Lightweight and durable aluminium frame
- Available frame colour: Silver/Red





# 2GOability Glide Pro

## The Compact Foldable Wheelchair

The 2GOability Glide Pro wheelchair is a high quality self-propelled, lightweight aluminium wheelchair. Weighing only 11.5 kg and featuring a folding frame, to allow for easy and convenient transportation. This stylish yet comfortable wheelchair is the definition of ideal.

- Folding backrest
- Height adjustable armrests
- Quick release fold in one mechanism engaged by pulling the seat upholstery
- Available frame colour: Silver



# Specifications

2GOability Pace II   2GOability Access   2GOability Glide Pro

Total Length	75 cm	90 cm	95 - 104 cm*
Total Width	62 cm	54 cm	67 cm
Total Height	101 cm	92 cm	95 cm
Seat Width	37 cm	45 cm	45 cm
Seat Depth	24 cm	43 cm	43 cm
Seat Height (From Ground)	55 - 65 cm*	50 cm	50 cm
Backrest Height		42 cm	44 cm
Total Weight	6.4 kg	9.5 kg	11.5 - 15 kg**
Maximum User Weight	100 kg (15 Stone)	100 kg (15 Stone)	113 kg (17 Stone)
Maximum Basket Weight	5 kg		
Adjustable Push Handles	Yes		
Push Handles Width	52 cm		
Transit		Yes	No
Self-propel		No	Yes
Backrest Foldable		Yes	Yes
Back Angle Adjustment		No	No
Attendant Brakes		Yes	No
Anti-tip Wheels		No	No
Armrests (Swing-away)		No	No
Armrests (Removable)		No	No
Armrests (Adjustable)		No	Yes
Extra Wide Comfort Pads		No	No
Height Adjustable Push Handles		No	No
Seat Upholstery (Removable)		No	Yes
Padded Seat		No	Yes
Legrests (Removable)		Yes	Yes
Calf Strap/Heel Strap/Calf Plate		Calf Strap	Heel Strap
Crash Tested		No	Yes
Available Colour Options	Silver	Silver/Red	Silver

**KEY:** \* = Adjustable or optional \*\* = Varies on prescription options



[www.vanosmedical.co.uk](http://www.vanosmedical.co.uk)

Van Os Medical UK Ltd | Excel House | Ashbrooke Park  
Sherburn-In-Elmet | North Yorkshire | LS25 6PJ

T : 01977 681 400 E : [sales@vanosmedical.com](mailto:sales@vanosmedical.com)