

Smart Ramp

- » The Smart Ramp has been designed for use with wheelchairs, powerchairs and scooters
- » Can be used to obtain safe access to motor vehicles, buildings and other inaccessible locations
- » A generous 76cm (30") width makes this ramp suitable for MPVs of up to 70cm (27.5") in width
- » When folded out they are 184cm (6ft) long
- » Lightweight aluminium construction
- » Folds in half for ease of storage and handling
- » Carries like a suitcase when not in use
- » Anti-slip high traction surface
- » Striking metallic colour finish
- » Weight capacity: 200kg (32st)


CODE: R008



Folding Ramp

- » The Folding Ramp provides a perfect and safe solution for disability and mobility access
- » Available in two sizes - 3ft and 6ft
- » Ideal for Wheelchairs, Powerchairs and Scooters
- » The two hinged tracks simply unfold and can be positioned accordingly for maximum safety and support with an angled lip at the top of the ramp for added security
- » Ideal for use by retail shops on the high street in accordance with the DDA (Disability Discrimination Act)
- » For temporary or permanent installation - two holes included for permanent installation
- » Features raised edges for added safety and peace of mind
- » These folding aluminium ramps are versatile, lightweight and easy to use
- » Constructed from lightweight aluminium and includes a central plastic handle for easy carrying and transportation
- » Two aluminium sections 35cm and 37cm wide to fold to a compact size
- » High quality user friendly alloy ramp
- » A wide anti-slip grip surface ensures safe use in most weather conditions
- » Durable design to support both wheelchair/mobility scooters and the occupant with a weight capacity of up to 272kg (43st)


CODE: R009



Channel Ramp

- » This 7ft Channel Ramp has been designed for use with manual wheelchairs, powered wheelchairs or scooters
- » It allows a wheelchair to travel up steps, over door thresholds or up into the back of vehicles
- » Easy to use and lightweight, this telescopic two-piece ramp extends to 2.1m (7ft) in length
- » to accommodate use in a variety of settings, then retracts again for ease of storage and transportation
- » A slip-resistant coating improves wheelchair grip
- » Folded dimensions: 100cm (39.4") x 26cm (10.2") x 6.5cm (2.6")
- » Ramp bag available for storage and transportation
- » Weight capacity: 272kg (43st)


CODES: R003 (Ramp)
RA003 (Bag)



Roll Up Ramp

- » This excellent 3ft Roll Up Ramp provides a quick and easy solution to obtain access to buildings and vehicles
- » It has been designed for scooters and wheelchairs and can be used as a portable or fixed ramp
- » Lightweight aluminium and rolls up to a compact size
- » Comes complete with a bag for easy storage and transportation
- » Weight capacity: 300kg (47st)


CODES: **RO04 (Ramp and Bag)**
RA004 (Bag)



Threshold Ramp

- » Fixed height design means the ramps can be simply placed over most doorways and can be easily removed when required
- » Lightweight with an excellent check plate grip surface
- » Portable - no fitting/ installation required
- » Aluminium finish prevents stains and corrosion
- » Comes complete with anti-slip pads on the underside of the ramp to offer protection to the floor surface whilst preventing slippage
- » Available in three sizes- 60cm (1.9ft), 80cm (2.6ft) and 100cm (3.2ft)
- » 70cm (27") wide - suitable for most scooters and wheelchairs
- » Narrow enough for most doorways
- » Maximum user weight 272kg (42st)

CODE: **RO10**



Ramp Dimensions

CODE	DESCRIPTION	LENGTH	WIDTH	LENGTH WHEN FOLDED	WIDTH WHEN FOLDED	WEIGHT	WEIGHT CAPACITY
RO08	Smart Ramp	6ft (1.8m)	76cm (30")	184cm (72")	38cm (15")	19.4kg (42lb)	200kg (32st)
RO09-3FT	Folding Ramp	3ft (0.9m)	72cm (28")	91cm (36")	32cm (12.5")	9kg (20lb)	272kg (43st)
RO09-6FT	Folding Ramp	6ft (1.8m)	72cm (28")	183cm (72")	32cm (12.5")	18kg (40lb)	272kg (43st)
RO03	Channel Ramp	7ft (2.1m)	53cm (21")	100cm (39.4")	26cm (10.2")	14.5kg (32lb)	272kg (43st)
RO04	Roll Up Ramp	3ft (0.9m)	76cm (30")	N/A	N/A	8kg (17.5lb)	300kg (47st)
RO10-600	Threshold Ramp - Small	1.9ft (0.6m)	59cm (23")	N/A	N/A	4kg (9lb)	272kg (42st)
RO10-800	Threshold Ramp - Medium	2.6ft (0.8m)	70cm (27")	N/A	N/A	5kg (11lb)	272kg (42st)
RO10-1000	Threshold Ramp - Large	3.2ft (1m)	70cm (27")	N/A	N/A	6kg (13lb)	272kg (42st)

Leaflet Code: LL230

For more information contact your local dealer:

All sizes and weights stated are nominal.
All errors and omissions excepted.
Due to the limitations of photography and the printing process colours shown may not be 100% accurate.

