

EXCEL G-LOGIC



www.vanosmedical.co.uk

The Excel G-Logic

Van Os Medical's greatest source of pride was the introduction of the markets first modern and inclusive lightweight aluminium manual wheelchair; the G-Lite Pro. Designed with aesthetics and function, working together we brought to you a chair that not only looked great, but was durable and simple to use.

Move forward 10 years and Van Os bring you the latest incarnation in the G-Series, The Excel G-Logic. Using the same proven chassis and frame, with enhancement of the wheelchair, by adding a new rear wheel inspired by the cycling industry, with a sleek aluminium hub, lightweight rim and push rim to give the chair a more sporty feel as well as look.


Mix and Match with a choice of four frame colours and four upholstery colours, its so easy to personalise your G-Logic and make it your own. Cleverly designed cushions with velcro strips also ensure you can add or remove cushions with no hassle even on the go.



Features

- Half-folding back
- Puncture proof tyres
- Height adjustable armrests
- Available in Self-Propel or Transit
- Interchangeable upholstery colours
- Optional elevating legrests available
- Swing-away and quick release legrests
- Anti-tip wheels as standard (Self-Propel)
- Quick release wheels for transportation (Self-propel)
- Attendant brakes as standard for both Transit and Self-propel

Available Frame Colours



Upholstery Colours



Easy to remove
Attendant brakes



Anti-tip wheels as
standard



Transit version
available

New and Improved

With the addition of a choice of four frame colours; Ice White, Silk Silver, Matt Black (Standard) or Trail Black* (UK Exclusive)* and four upholstery colours, the Excel G-Logic is now even more customisable than ever with a possibility of 16 variations.

Structural redesigns include a new side frame that now includes a new 4 degree seat angle, providing the most comfortable and ergonomic posture. The backrest angle of the G-Logic has also been adjusted to give the user a more natural seating position, which now aligns to the curvature of the spine to give superb lumbar support.

The Excel G-Logic is also Crash Test approved.

G-Logic 



Foldable

Easily foldable for transportation, the G-Logic can fit into the boot of almost any car. This wheelchair is the perfect travel companion for frequent travellers with an active lifestyle.

Featuring a quick release mechanism engaged by pulling up the seat upholstery, folding the chair up in one easy simple motion; convenience couldn't be easier.

- Half-folding back
- Removable seat cushion
- Swing-away and quick release legrests
- Quick release wheels for transportation
- Foldable footplates for more compact storage



G-Logic Specifications

Seat Width	40, 45 , 50 cm
Seat Depth	43 cm
Seat Height	50 cm
Backrest Height	44 cm
Total Length	95 - 104 cm*
Total Width	62, 67, 72 cm
Total Height	95 cm
Total Weight	12 - 15 kg**
Maximum User Weight	113 kg (40 & 45 cm) 135 kg (50 cm)
Transit	Yes
Self-propel	Yes
Backrest Foldable	Yes
Back Angle Adjustment	No
Attendant Brakes	Yes
Anti-tip Wheels	Yes*
Armrests (Swing-away)	No
Armrests (Removable)	No
Armrests (Adjustable)	Yes
Extra Wide Comfort Pads	No
Height Adjustable Push Handles	No
Seat Upholstery (Removable)	Yes
Padded Seat	Yes
Legrests (Removable)	Yes
Calf Strap/Heel Strap/Calf Plate	Heel Strap
Crash Tested	Yes
Available Colour Options	White, Silver, Matt Black, Trail Black***

KEY:

Bold = Standard option

* = Adjustable or optional

** = Varies on prescription options

*** = UK Exclusive





Van Os Medical UK Ltd | Excel House | Ashbrooke Park
Sherburn-In-Elmet | North Yorkshire | LS25 6PJ

T : 01977 681 400 E : sales@vanosmedical.com

Follow us



www.vanosmedical.co.uk