

# TITAN AXS

“Manoeuvrable and intelligent addition to the Titan family”



Designed for optimal drive traction on gradients

Red & blue interchangeable colour options

CE

The Titan AXS is a manoeuvrable powerchair which can be used both indoors and outdoors. Compact and agile, the Titan AXS moves effortlessly around tight corners, it's ideal in areas where space is limited.

[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

**drive**  
DeVilbiss  
HEALTHCARE

The AXS is the latest mid-wheel addition to the Titan family. It has been configured to provide maximum performance, both indoors and outdoors, thanks to its drive wheel position and flexible front castor wheels.

## Features

- › Maximum speed of 4mph\*
- › Maximum range of 19 miles on a full battery charge\*
- › Maximum weight capacity 136kg (21st)
- › Mid-wheel drive for optimum manoeuvrability
- › Six wheels provide added stability when using the powerchair
- › Height adjustable seat with high back and headrest for added comfort
- › Swivel seat facilitates easy transfers on and off the chair
- › Freewheel Mode to allow use of the powerchair without switching the motor on
- › Height adjustable and removable headrest to suit user requirements
- › Flip-up armrests allow improved access to tables, desks etc. and allowing to facilitate easy transfers
- › Dynamic LiNX programmable controller provides an easy and enjoyable driving experience
- › Moulded PU armrests that are width, height and angle adjustable to suit comfort preferences
- › 25cm (10") middle wheels and 15cm (6") front and rear castors with solid puncture proof tyres provide durability and a low maintenance
- › Angle and height adjustable flip-up footplate to facilitate transfers
- › Powerful 320W motors
- › Rear anti-tip wheels as standard
- › Seatbelt included for additional safety and security

\*Subject to conditions of use


Colour Options:

Red

Blue



New Intelligent Dynamics LiNX controller



Well padded, reclining & height adjustable seat

## Key Features

- Maximum range of 19 miles on a full battery charge\*
- Easily programmable controller with Bluetooth facility
- Mid wheel drive for optimum manoeuvrability
- Powerful 320W motors



Configured for manoeuvrability

## Specifications:

Max Weight Capacity	135kg (21st)
Max Speed	4mph (6.4km)
Estimated Range	19m (30.5km)
Max Climbing Angle	6°
Turning Radius	51cm (20")
Ground Clearance	6cm (2.5")
Powerchair (L x W x H)	100cm x 59cm x 112cm (39.25" x 23.5" x 44")
Seat (L x W x H)	46cm x 45cm x 46cm (18" x 18" x 18")
Armrest to Seat	18cm
Backrest Height, excl. Headrest	42cm
Height (Folded Backrest)	77cm
Castors	15cm x 5cm (6" x 2")
Mid Drive Wheels	25cm x 8cm (10" x 3")
Motor	320W x 4600 RPM
Controller	Dynamic LiNX 40A
Batteries	12V x 33Ah
Battery Charger	4A
Base Weight	40kg (88lb)
Battery Weight (Pair)	20kg (44lb)
Seat Weight	18kg (40lb)
Total Weight	78kg (172lb)

Leaflet Code: LL335

If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.