




COMFORT EDITION WALKING STICK

This innovative new walking aid retains all the popular features of the original HurryCane but also offers some additional benefits that make it this HurryCane Walking Stick..... The Comfort Edition!!



New to the HurryCane Comfort Edition

- » New pastel colours including blue, green and pink, as well as black
- » New and improved moulded handle to fit the hand more comfortably
- » New design three tip base helps the stick stand independently

Features:

The HurryCane gives three points of contact for superior balance and stability. It also pivots to help the user maintain balance on inclines, uneven ground and shifting terrain.

- » Free standing design
- » 360° pivoting head
- » Folds easily for storage
- » Height adjustable from 77.5cm to 95cm
- » 8 convenient height adjustments
- » Anodized aluminium shaft
- » Wrist strap & travel bag included
- » Maximum user weight 159kg

Specifications:

Length 77.5-95cm (30.5-37.4")	Product Weight 430g (0.9lb)	User Weight Capacity 159kg (25st)
----------------------------------	--------------------------------	--------------------------------------

Codes:

Code	Description
HCANE2-BLK	HurryCane Comfort Edition - Black
HCANE2-PBL	HurryCane Comfort Edition – Pastel Blue
HCANE2-PGR	HurryCane Comfort Edition – Pastel Green
HCANE2-PPK	HurryCane Comfort Edition – Pastel Pink

drive

DeVilbiss
HEALTHCARE

www.drivedevilbiss.co.uk