

AIDAPT BATHROOMS LTD.Lancots Lane
Sutton Oak
St. Helens
WA9 3EX

The Consort Deluxe Self Propelled Steel Wheeled Chair


The Consort Deluxe Self Propelled Steel Wheeled Chair from Aidapt. Made from high grade durable steel with a double cross bar for strength. The Consort offers ease of movement for the user and quick-folding mechanism for trouble-free transportation. Features include quick release rear wheels, flip up armrests, 20" wide padded seat for comfort, hand brakes for the carer, anti tilt wheels at the rear, PU armrests, foldable backrest and detachable footrests. Maximum user weight limit 180kg (28 stone).

Features:

- Flip up armrests
- 20inch seat width (510mm)
- Detachable footrests
- Foldable backrest
- Hand brakes
- Quick release rear wheels

PRODUCT SPECIFICATION

Length (cm):	1060
Height (mm):	955
Width (mm):	710
Product Dimensions (mm):	955x710x1060
Net weight (kg):	20
Packaged dimensions (HxWxL mm):	800x360x760
Display Footprint (mm):	955x710x1060
Transport dimensions (mm):	800x360x760
Gross Weight (Packaged) (kg):	21.5
Packaging:	Plain Carton
Item priced per:	1
Warranty Period:	12 Months
Backrest Height (mm):	410
Colour:	Silver
Configuration:	Self Propelled
Front Wheel (mm):	200
Maximum User Weight (kg):	180
Maximum User Weight (stones):	28
Overall Height (mm):	955
Overall Width (mm):	710
Rear Wheel (mm):	610
Seat Depth (mm):	415
Seat Dimensions (mm):	510x415
Seat Height (mm):	520
Seat Width (mm):	510
Size:	20"
Width (Extended) (mm):	710
Width (Folded) (mm):	320



Product Code: VA166SILVER

Barcode: 5021196800859

Registered Office:
AIDAPT BATHROOMS LTD.
 LANCOTS LANE, SUTTON OAK, ST. HELENS, WA9 3EX