

A Series Rollators

Description	Product Code
A Series Rollator - Red	P11255
A Series Rollator - Blue	P11267
A Series Rollator - Silver	P11279
A Series Petite Rollator Grey	S15232
A Series Tall Rollator Grey	S23472

P11267

S23472



Features and Benefits

- Different colour and height options to suit different users
- Backrest and comfortable seat to allow use as a chair
- Simple folding mechanism for easy storage and transport
- Large diameter wheels with high-grip tyres and lockable brakes for safety
- Under-seat storage bag for convenience

A Series Rollators

Product Description

A choice of lightweight rollators to help users who have difficulty or lack confidence when walking. The regular options are suitable for medium height users and are available in three colours - red, blue and silver. The petite and tall options are suited to smaller and taller users and are available in grey. All options also have adjustable handle height for more tailored size changes to suit the needs of the user.

Brake levers on each handle control the brakes on the two rear wheels. They also include a locking mechanism for added safety. The large diameter high-grip wheels allow the rollator to roll along most surfaces while keeping control. A cushioned backrest and padded seat allow use of the rollator as a chair. The simple folding mechanism allows the user to easily collapse and store the rollator, useful when loading into a car boot or storage cupboard.

Each rollator includes a removable under-seat storage bag to hold smaller items, keeping hands free for safety.

Product Technical Details

- A Series Petite (S15232) for users 5'0" to 5'6".
- A Series Tall (S23472) for users 5'6" to 6'3".
- Painted aluminium frame with plastic parts.
- Lockable brakes to rear wheels.
- Removable storage bag, capacity 5 kg.
- Maximum User Weight: 135 kg (21 st).
- Product weights: P11255, P11267, P11279 and S23472: 6.5 kg. S15232: 6.3 kg.

All Models



A Series Rollators

P11255, P11267 and P11279



S15232





S23472



Standard Features

• NRS guarantee - refer to the Terms and Conditions on our website

NRS Healthcare. Sherwood House, Cartwright Way, Forest Business Park, Bardon Hill, Coalville, Leicestershire, LE67 1UB, UK Tel: +44 (0)345 121 8111 Email: <u>customerservice@nrshealthcare.co.uk</u> Web: <u>www.healthcarepro.co.uk</u>