

drive



aerocarbon

POWERCHAIR

Weighing just

12kg

Carbon Fibre
Technology

Warranty

5 Years

drivedevilbiss.co.uk

AeroCarbon POWERCHAIR

Carbon fibre technology delivers exceptional strength without the weight. As part of the Aero Carbon family, this sleek, modern powerchair is designed for effortless manoeuvrability, advanced comfort, and intuitive control with optional attendant control.

Specifications

| | |
|------------------------------|---------------------------|
| Length | 91cm (35.8") |
| Width | 61cm (24") |
| Height | 89cm (35") |
| Seat Depth | 42cm (16.5") |
| Seat Width | 44cm (17.3") |
| Seat Height | 44cm (17.3") |
| Backrest Width | 42cm (16.5") |
| Folded Size | 73x31x76cm (LxWxH) |
| Ground Clearance | 3.2cm (1.25") |
| Battery Supplied | 24V 10Ah Lithium-ion |
| Motor Type | 24V / 150W (DC Brushless) |
| Range* | 9.5 miles |
| Maximum Speed* | 4mph |
| Carry Weight (excl. battery) | 12kg |
| Total Weight (incl. battery) | 13.5kg (29.8lb) |
| Maximum User Weight | 136kg (21.4 stone) |

AeroCarbon Powerchair:
ULCFPCHAIR

Attendant Controller & Fixings:
12-EC8003Z-00

WARRANTY:
Frame - 5 years.
Electrical Components - 3 years.
Battery - 1 year.
Excludes wear & tear, & misuse.



drive

£

Contact us for more information

Our lightest
powerchair
weighs in at
just 12kg



Perfect for
transport, with
easy folding
and ideal for
car boots



Premium quality,
lightweight
carbon fibre
frame with a 5
year warranty



Air travel safe,
lightweight
lithium battery
with off board
charging



Optional
easy-to-fit
attendant
control
available

