

4 Wheel Folding Forearm Rollator

Description	Product Code
4 Wheel Folding Forearm Rollator	S20756



Features and Benefits

- Foam forearm pads for comfort
- Lightweight, foldable aluminium frame for easy transport
- Adjustable handle height and distance to handles to suit different users
- Double brakes with additional parking lock
- Includes a seat and backrest to allow the user to rest during use
- Walking stick holder and Storage Bag for convenience
- Large diameter high-grip wheels to cope with a range of floor surfaces

4 Wheel Folding Forearm Rollator

Product Description

A lightweight rollator to help users who have difficulty or lack confidence when walking. This rollator also has forearm supports to help users that may have limited grip strength or upper body dexterity. A foam-covered backrest and a seat allow use of the rollator as a chair. The backrest fits using button-clips and is key to adding stability when assembled, but easily disconnects to allow folding the rollator for storage. The lockable brakes act on the rear wheels. Simple button-clips and handwheels allow the user to adjust the height of the forearm supports and distance to handles. The rollator has a walking stick clip and holder to one side and a detachable Storage Bag to the front to help the user to easily transport lightweight objects. The large diameter high-grip wheels allow the rollator to roll along most surfaces while keeping control.

Note: We recommend that a healthcare professional assesses a user before they buy and use this product.

Product Technical Details

- Mostly painted aluminium construction. Cushioned forearm supports and foam cover to backrest.
- Detachable Storage Bag for small items.
- Lockable brakes on rear wheels.
- Product Weight: 9.2 kg (1.5 st).
- Maximum User Weight: 136 kg (21 st).



Standard Features

- NRS guarantee - refer to the Terms and Conditions on our website

NRS Healthcare. Sherwood House,
Cartwright Way, Forest Business Park,
Bardon Hill, Coalville,
Leicestershire, LE67 1UB, UK

Tel: +44 (0)345 121 8111
Email: customerservice@nrshealthcare.co.uk
Web: www.healthcarepro.co.uk