

**drive**



# aerowalk

## carbon

---

**4.6kg\***  
Weighing just

---

---

Lightweight material  
**Carbon Fibre**

---

---

**Travel Ready**  
Compact Folding

---

[drivedevilbiss.co.uk](http://drivedevilbiss.co.uk)

# AeroWalk Carbon

Part of the Aero Carbon Family, the AeroWalk Carbon Rollator is the perfect blend of cutting-edge design, strength, and effortless mobility. Crafted from premium carbon fibre, it offers exceptional durability while remaining ultra-lightweight for easy handling and transport.

## Specifications

Length	68cm (26.8")
Width	60cm (23.6")
Folded Length	68cm (26.8")
Folded Width	21cm (8.3")
Folded Height	84cm (33.1")
Seat Depth	22cm (8.6")
Seat Width	45cm (17.7")
Seat Height	57cm (22.4")
Handle Height	80 – 95cm (31.5 – 37.4")
Front Wheel	9"
Rear Wheel	8"
Width Between Hand Grips	48cm (18.9")
Product Weight (excl. accessories)	4.6kg (10.1lb)
Product Weight	5.1kg (11.2lb)
Maximum User Weight	150kg (23.6 stone)

AeroWalk Carbon Fire Red:  
720200030

AeroWalk Carbon Butter Cream:  
720200040

Fire  
Red

Butter  
Cream

**drive**

£

Contact us for more information

Our lightest  
rollator  
weighs in at  
just 4.6kg



Perfect for  
travel with  
compact  
folding



Premium quality  
carbon fibre  
material provides  
strength & durability  
while remaining  
lightweight



Comfortable  
seat and  
ergonomic  
handles



All round  
reflectors for  
enhanced  
safety

