

drive



EasyFold

FOLDING ELECTRIC WHEELCHAIR



- **Easily folds down for travel & transport -
Folded measurements 73 x 36 x 77 cm (28.7 x 14.2 x 30.3")**
- **Travels up to 12 miles on a full charge at speeds of up to 4mph**
- **Off board charging facility allows you to charge the batteries independently**
- **Comes complete with a removable cushion for extra comfort**

The Easy-Fold electric wheelchair can be easily folded & unfolded for travel or transport with a manageable carrying weight of 36kgs, the battery & footrests can be removed as well as a folding backrest for a more compact size which can be easily stored at home, in your car boot or loaded onto public transport. (73 x 36 x 77 cm (28.7 x 14.2 x 30.3"))

The Easy-Fold mobility powerchair has adjustable and flip up armrests to allow for an easier transfer to and from your powerchair as well as swing away footrests to help you gain a better seating position at the table, making it possible to be closer to your family & friends.

- Length: 106cm (42 ")
- Width: 66cm (26")
- Height: 90cm (36")
- Seat Depth: 42cm (16.5")
- Seat Width: 46cm (18")
- Seat Height: 54cm (21")
- Max User Weight: 136kg (21st)
- Carry Weight: 36kgs (79lbs)
- Battery Capacity: 2 x 34AH
- Battery Weight: 22kg (48lb)
- Range: 12 miles (19km)